


Title	CD34 <sup>+</sup> cell reprogramming using episomal vectors
Date Submitted	May 5, 2012
Submitted by -	Olga Momcilovic, Zeng Lab <a href="mailto:omomcilovic@buckinstitute.org">omomcilovic@buckinstitute.org</a>
Adapted from -	
Contributors -	Olga Momcilovic
Affiliation(s) -	Zeng Lab Buck Institute

## ❖ Introduction:



iPSC colony derived from CD34<sup>+</sup> cord blood cells by reprogramming with non-integrating plasmids.

## ❖ Protocol:

1. Prime CD34<sup>+</sup> cells
  - a. Thaw CD34<sup>+</sup> cells using Lonza's protocol and culture for 4-5 days ([http://www.lonzabio.com/uploads/tx\\_mwaxmarketingmaterial/LonzaManualsProductInstructionsProcedureforThawingPoieticsCells.pdf](http://www.lonzabio.com/uploads/tx_mwaxmarketingmaterial/LonzaManualsProductInstructionsProcedureforThawingPoieticsCells.pdf))
2. Day 1: nucleoporate 1x10<sup>6</sup> hCD34<sup>+</sup> cells with single (up to 10 μg), or combination of plasmids (8 μg C5 + 2 μg Tg) by Amaxa using program U-008

3. Days 1 and 2: culture nucleoporated cells in one well of a 12-well plate in the CD34<sup>+</sup> medium with cytokines
4. Day 3: transfer nucleoporated cells to 3 wells of MEF coated 12-well plate and culture in MEF medium for one day
  - a. Once cells are seeded into wells, spin plates at 100xg for 30 min to help cells attach to MEF coated wells
5. Day 4: replace MEF medium with hESC medium (supplemented with 10 ng/ml FGF2)
  - a. OPTIONAL: collect MEF medium and spin it down at 100xg for 5 min; aspirate medium and resuspend cell pellet in hESC medium with 10 ng/ml FGF2 – some CD34<sup>+</sup> cells may not attach during the first day, so save them and replat them
  - b. Change hESC medium every other day for total of 6 days
  - c. OPTIONAL: add valporic acid (0.5 mM) or Na-butyrate (0.25 mM)
6. Switch to MEF-CM with 10 ng/ml FGF2 one week after transfer onto MEF coated wells
7. Two weeks after nucleoporation, perform TRA1-60 staining on live cells to identify most likely iPSC clones
  - a. With cord blood CD34<sup>+</sup> cells expect to see colonies appearing 7-11 days post-nucleoporation
  - b. With adult bone marrow and peripheral blood CD34<sup>+</sup> cells colonies start appearing 11-14 days post-nucleoporation
8. Manually dissect each TRA1-60 positive colony and transfer to a separate well of a 12-well plate: each colony becomes a clone
  - a. OPTIONAL: add 10  $\mu$ M ROCK inhibitor and/or hESC cloning and recovery supplement to improve survival and attachment of dissected colonies
9. For the first 2 – 3 passages keep clones in 12-well plates, then expand to 35 mm dishes
10. Manually passage clones for the first 6 – 10 passages, then switch to 1 mg/ml collagenase (depending on whether clones remain undifferentiated when enzymatically passaged). In instances when less than 10% of colonies are differentiated, remove differentiated cells manually and proceed to enzymatic passage; if more than 10% colonies are differentiated, continue with manual passaging
11. Gradually reduce FGF2 concentration in MEF-CM to 4 ng/ml and switch to hESC medium by mixing MEF-CM and hESC medium in order to adopt iPSC clones to hESC medium with 4 ng/ml of FGF2.

❖ **Materials:**

Product	Company	Catalogue number
MEF, mitomycin C treated	Millipore	PMEF-N
DMEM, high glucose	Gibco	11995
FBS		

KNOCKOUT™ DMEM/F12	Gibco	12660
NEAA	Gibco	11140
Anti-Anti	Gibco	15240
KNOCKOUT™ Serum Replacer	Gibco	10828
2-mercaptoethanol	Gibco	21985
GlutaMAX™-1	Gibco	35050
CD34+ cells	Lonza	2C-101
HPGM™	Lonza	PT-3926
DNase I	Sigma	D4513
SCF	Peprtech	AF-300-07
TPO	Peprtech	AF-300-18
FL	Peprtech	AF-300-19
FGF2	Stemgent	03-0002
Nucleofector kit for CD34+ cells	Lonza	VPA-1003
ROCK inhibitor, Y27632	Stemgent	04-0012
hESC cloning and recovery supplement	Stemgent	01-0014-500
Na-butyrate	Stemgent	04-0005
Valporic acid	Stemgent	04-0007
TRA1-60 antibody	eBioscience	13-8863-83
C5 - EBNA1 carrying OCT4, SOX2, KLF4, LIN28, cMYC	Addgene	<a href="http://www.addgene.org/28213/">http://www.addgene.org/28213/</a>
TG - EBNA1 carrying SV-40 Large T antigen	Addgene	<a href="http://www.addgene.org/28220/">http://www.addgene.org/28220/</a>

MEF medium  
90% DMEM  
10% FBS  
1% Anti Anti

hESC/hiPSC medium  
KNOCKOUT™ DMEM/F12  
20% KNOCKOUT™ Serum Replacer  
1% GlutaMAX™-1  
1% NEAA  
1% Anti/Anti  
4 - 10 ng/ml FGF2  
0.1 mM 2-mercaptoethanol

CD34+ cell medium (recommended by Lonza)  
HPGM™ Hematopoietic Progenitor Growth  
Medium supplemented with the following  
concentrations of cytokines:  
FL - 50 ng/ml

TPO – 50 ng/ml

SCF – 25 ng/ml

All cytokines are from Peprotech and are diluted in trehalose at concentration of 100 ng/ $\mu$ l.

Abbreviations

MEF = mouse embryonic fibroblasts

FBS = fetal bovine serum

NEAA = non-essential amino acids

FL = Flt3 ligand

SCF = stem cell factor

TPO = Thrombopoietin

MEF-CM = hESC medium conditioned for 24 hrs on MEF

❖ Troubleshooting:

❖ **References:**