

Carbonate Extraction of Nuclear Envelopes (Rout Lab, 1999)

Purpose: to extract proteins from the nuclear envelope using carbonate

1. Mix 25 μ l of a nuclear envelope fraction with 225 μ l BT/Mg/DMSO in TLA-45 tubes. Less NE can be used if necessary and brought up in a total volume of 250 μ l with BT/Mg/DMSO.
2. Spin tubes in TLA-45 rotor for 40K, 40 min., 4°C.
3. Aspirate supernatant and resuspend remaining pellet with 250 μ l of 0.1M Na-Carbonate extraction buffer (this is usually a 1:10 dilution of a 1M Na-Carbonate buffer stock with water). Vortexing and pipetting up and down will be sufficient to resuspend pellet. Do NOT sonicate.
4. Incubate tubes on ice for 1 hour.
5. Spin tubes in TLA-45 rotor for 40K, 40 min., 4°C.
6. Transfer supernatant to a fresh tube labelled "S". Wash pellet with 250 μ l Na-Carbonate extraction buffer.
7. Spin tubes in TLA-45 rotor for 40K, 40 min., 4°C.
8. Transfer supernatant to "S" tube. Add 500 μ l water to "S" tube and 750 μ l water to tube containing pellet.
9. Add 100 μ l 0.3% NaDOC and 100 μ l 72%TCA to all tubes to precipitate proteins. Mix by inversion and incubate on ice for 1 hour.
10. Spin tubes in microcentrifuge at 14K for 20 min., 4°C.
11. Aspirate supernatant and resuspend pellet in 0.5 mL 90% acetone. Incubate in freezer for 3 hours to overnight. (Overnight is preferred)
12. Spin tubes in microcentrifuge at 14K for 20 min., 4°C.
13. Aspirate supernatant and resuspend pellet in 5 μ l sol. A and 5 μ l sol. B. Heat at 65°C for 10 minutes and run on protein gel.