

Day 1

DATE: 01-20-09

1 ChIP of *** cells: Preparation of beads, binding primary Ab, followed by chromatin

PURPOSE: Isolation of modified histones from chromatin

MATERIALS & METHODS:

Antibodies (5ug/IP):

Chromatin (use 0.5-0.35mg of chromatin in 1mL):

10% DOC (sodium deoxycholate –Sigma-Aldrich)

5mg/mL BSA (0.125g BSA in 25mL 1X PBS) make immediately before use, put on ice

Sheep-anti-rabbit beads (antibody coupled magnetic beads - Dynal/Invitrogen)/50ul per IP

NOTE: larger volume of beads can be washed at one time for convenience

#	Epitope	Cat#	Lot#	Species	[Ab]	Ab amt.	Vol/IP	Vol PBS
T01	H3K4me3	ab8580	500368	rabbit	1mg/ml	5ug	5ul	995
T02	H3K4me3	ab8580	500368	rabbit	1mg/ml	5ug	5ul	995
T03	H3K4me1	ab8895	133163	rabbit	.5mg/ml	5ug	10ul	990
T04	H3K27me3	07-449	Dam1514011	Rabbit	1 mg/ml	5ug	5ul	995
T05	H3K36me3	9050	Dam475787	Rabbit	1mg/ml	5ug	5ul	995
T06	CTCF	9Mix	1/23/09	mouse			25ul	975

1. Use magnetic strip. Remove supernatant with a pipette/vacuum and resuspend in 1mL PBS containing 5mg/ml BSA to each tube (make immediately before use from Sigma BSA powder, cat. A-3350).

2. Wash 3 times. Resuspend beads in 9XXuL PBS+BSA per tube (see chart).

3. With tube against the magnet, add antibody for a total volume of 1000ul. Incubate for **6-8hrs** on a rotating platform at 4°C. **NOTE:** if antibody is being replicated, 100ul of beads and 10ug of antibody can be combined in 1000ul for binding, split ab/beads in half before adding chromatin below.

Make 10% DOC (0.1g DOC in 1mL dH2O) immediately before use

Make 5mg/mL BSA (0.125g BSA in 25mL PBS) immediately before use

Make 50X Complete (1 tablet in 1.0mL dH2O)

4. Remove supernatant with a pipette/aspirate and resuspend in 1mL PBS+BSA (on ice)

5. Wash 3 times. (let sit on magnet)

6. Resuspend antibody beads in 100uL 1X PBS+BSA

7. Set up IP reactions with chromatin extract.

#	Chromatin	Description
T01	XX (0.4mg/ml)	H3K4me3
T02	XX (0.4mg/ml)	H3K4me3
T03	XX(0.4mg/ml)	H3K4me1
T04	XX (0.4mg/ml)	H3K27me3
T05	XX(0.4mg/ml)	H3K36me3
T06	XX(0.4mg/ml)	CTCF

Add following additional reagents to bring the volume to 1.3 mL with following concentration:

Stock Solution	Final Concentration	Volume Per Tube(ul)	7 X (ul)
10% Triton-X	1%	130	910
10% DOC	0.10%	13	91
50X Complete	1X	26	182
1X TE	1X	131	917

300ul per tube

8. Add 1.3mL chromatin mix to magnetic beads and incubate at 4°C overnight in a rotating platform.

Day2

DATE: 01-22-09

2 ChIP of *** cells: Washing beads, Elution, and reversal of crosslink

Make 10% DOC (1g DOC in 10mL dH₂O/.35g in 3.5ml/.7g in 7.0ml) immediately before use.

Make RIPA buffer immediately before use. Add the stock solutions in the order listed.

RIPA buffer:

Components	STOCK	50ml	100ml
dH ₂ O		34.8ml	69.6ml
50mM Hepes, pH 8.0	1M	2.5ml	5ml
1% NP-40	10%	5ml	10ml
0.7% DOC	10%	3.5ml	7ml
0.5M LiCl	8M	3.125ml	6.25ml
Complete	50X	1ml	2ml

Washing of Beads

1. Use a magnet MPC-E to precipitate the beads (optional: save the 1st supernatant). Keep tubes on ice. Wash **5 times** with 1mL RIPA buffer. Remove buffer by aspiration.

2. Wash once with 1mL TE.

3. After removing the TE by aspiration, quickly spin the tubes or place on magnet and remove remaining liquid with a pipet.

Elution from beads and reversal of cross-link

Elution buffer: (doesn't need to be fresh)

Components	STOCK	50mL
10mM Tris, pH 8.0	1M	500uL
1mM EDTA	0.5M	100uL
1% SDS	10%	5mL
dH ₂ O		44.35mL

1. Add 150uL of elution buffer, vortex briefly to resuspend the beads and incubate at 65°C for 10 minutes. Vortex briefly every 2 minutes during incubation, or put on thermomixer for 20min (**TURN ON AHEAD OF TIME**).

2. Spin for 30 seconds at maximum speed or place on magnet and transfer liquid to a new tube.

3. Reverse crosslink at 65°C O/N in incubator.

4. Also reverse crosslink 30uL input chromatin in 120uL elution buffer.

INPUT: XX

Day 3

DATE: 01-23-09

3 ChIP: Isolate - DNA precipitation

1. Add 250ul TE to each sample
2. Add 8ul of 10mg/ml RNase A (final conc. = 0.2mg/ml)
3. Incubate at 37°C for 1hr.
4. Add 8ul of 20mg/ml Proteinase K (final conc. = 0.4mg/ml)
5. Incubate for 1hr at 55°C.
6. Prepare Phase Lock tube per IP, spin at 14K rpm for 1min
7. Add 400ul Phenol:Chloroform: Isoamyl alcohol to each tube.
8. Add sample to Phase Lock tube and shake vigorously, do not vortex
9. Spin 4 min at 14K rpm. **Note:** if aqueous phase is cloudy, extract again.
10. Transfer aqueous layer to new tube
11. Add 16ul of 5M NaCl (final conc. = 200mM), plus 30ug of glycogen (1.5ul of 20mg/ml)
12. Add 920uL cold 100% EtOH, vortex briefly.
13. Incubate at -80°C for 20-30min. Spin at 14K rpm for 15 minutes at 4°C.
14. Wash pellet with 1mL cold 70% EtOH, vortex, spin 5min. at 4°C at 14K rpm.
15. Resuspend thoroughly in 34ul* of 10mM Tris (EB from qiagen kits). OD input DNA.
 - A) Proceed to Library Prep immediately for ChIP-Seq
 - B) Store at -80°C.